


Holding Pipettes

Quality – Made in Germany



Holding pipette for in-vitro fertilization and manipulation of the ovum

Features:

- Each batch is MEA (Mouse Embryo Assay) and LAL (Endotoxin Assay) tested
- Ideal in handling and safe for transportation
- With expedient size and angle
 - Lumen diameter set for stable immobilisation of cells
 - Optimal degree values
- Packaging
 - Individually sterile packed
 - Easy handling during pipette removal
 - Special brackets protect pipette from damage during transportation
 - Colour coding simplifies identification
- Validated electron beam sterilization
- Class Is according to Directive 93/42/EEC

Designed for precision work

Available in the following angles:

Holding pipette, angle 35° REF 14697

Holding pipette, angle 30° REF 14698

Technical data

Specification:

- Inner diameter $20 \mu\text{m} \pm 5 \mu\text{m}$
- Outer diameter $100 \mu\text{m} \pm 10 \mu\text{m}$

Holder:

- Color code GREEN
- Material POM

