



**Comments for pT-REx-DEST30**  
**7544 nucleotides**

CMV promoter: bases 1-588

TATA box: bases 502-508

Tetracycline operator (2X TetO<sub>2</sub>) sequence: bases 518-557

attR1 recombination site: bases 699-823

Chloramphenicol resistance gene: bases 932-1591

ccdB gene: bases 1933-2238

attR2 recombination site: bases 2279-2403

T7 promoter: bases 2464-2483 (complementary strand)

SV40 polyadenylation region: bases 2915-3045

f1 origin: bases 3175-3603

SV40 early promoter and origin: bases 3790-4098

Neomycin resistance ORF: bases 4157-4951

Polyadenylation region: bases 5015-5063

Ampicillin (*bla*) resistance ORF: bases 5474-6334

pUC origin: bases 6479-7152