

BIOL[®]GIX



66-1501 PETRI DISHES

EO Sterilization

PETRI DISHES



60mm (66-1560),
90mm (66-1501,
66-1701), or
150mm (66-1515)
in diameter

INOCULATING LOOPS & NEEDLES



1. The opposite end of each loop can be used as a needle
2. Gamma radiation sterilized, non-pyrogenic
3. Loops: 66-0001 (1 μ l), 65-0010 (10 μ l), needles: 65-0002

CELL SPREADERS



1. Smooth surface; L-shape
2. Gamma radiation sterilized, non-pyrogenic
3. Individual packaging (65-1001) or zip-seal packaging(65-1010)

www.BiologixGroup.com

inquiry@biologixgroup.com