

QuantiGene ViewRNA ISH Cell 740 Module

Cat. No. QVC0200

Product Overview

The QuantiGene[®] ViewRNA ISH Cell 740 Module is designed to work in conjunction with both the QuantiGene ViewRNA ISH Cell Assay Kit and QuantiGene ViewRNA miRNA ISH Cell Assay Kit and contains the signal amplification and detection components for visualization of target mRNA in the 740 channel using adherent or suspension cells. The QuantiGene ViewRNA ISH Cell 740 Module is supplied with enough reagents to process 24 assays when using the 24-well chamber plate assay format and 144 assays when using 96-well optical-bottom plates. The kit is configured for processing a minimum of 6-wells per experiment when using the 24-well plate assay format. If fewer assays are processed per experiment, reagent shortages may occur.

Intended Use

QuantiGene ViewRNA ISH Cell 740 Module is designed for use with both QuantiGene ViewRNA ISH Cell Assay Kit and QuantiGene ViewRNA miRNA ISH Cell Assay Kit and QuantiGene ViewRNA TYPE 10 Probe Sets for visualization of target mRNAs within individual cells.

Contents and Storage

The QuantiGene ViewRNA ISH Cell 740 Module contains the following components. Reagents have a minimum shelf life of 6 months from date of receipt when stored as recommended.

Table 1 QuantiGene ViewRNA ISH Cell 740 Module Contents and Storage:

Component	Description	Quantity	Storage
PreAmp10-740	DNA in aqueous buffered solution	420 µL	-20 °C
Amp10-740	DNA in aqueous buffered solution	420 µL	-20 °C
LP10-740	ALEXA 750-conjugated oligonucleotide in aqueous buffered solution	420 µL	-20 °C

Safety Warnings and Precautions

All chemicals should be considered potentially hazardous. We recommend that this product and its components be handled by those trained in laboratory techniques and be used according to the principles of good laboratory practice.

US Headquarters

Affymetrix, Inc.
3420 Central Expressway
Santa Clara, CA 95051
Tel: 1-888-362-2447
Direct: 1-408-731-5000
Fax: 1-408-731-5380
Email: info@affymetrix.com
Email: orders_fremont@affymetrix.com
Email: pqbhelp@affymetrix.com

European Headquarters

Panomics Srl
Via Sardegna 1
20060 Vignate-Milano (Italy)
Tel: +39-02-95-360-250
Fax: +39-02-95-360-992
Email: info_europe@affymetrix.com
Email: order_europe@affymetrix.com
Email: techsupport_europe@affymetrix.com

Asia Pacific Headquarters

Affymetrix Singapore Pte Ltd.
No.7 Gul Circle #2M-01/08
Keppel Logistics Building
Singapore 629563
Tel: +65-6395-7200
Fax: +65-6395-7300
Email: info_asia@affymetrix.com
Email: order_asia@affymetrix.com
Email: techsupport_asia@affymetrix.com

www.affymetrix.com/panomics

For research use only. Not for use in diagnostic procedures.