

BIOLIGIX®

PCR PRODUCTS-0.2ML TUBES

INDIVIDUAL PCR TUBES



Features:

- Made of medical grade polypropylene.
- Ultra thin-wall tubes provide faster heating and cooling.
- High degree of transparency allows easy observation of samples.
- Domed caps ensure efficient heat transfer, while flat caps are ideal for Real-Time PCR(qPCR).
- Unique design prevents opening during cycling and incubation.
- DNase, RNase, and endotoxin free.
- Fit most standard thermocyclers.
- Can withstand centrifugation of 6,000 RCF.




8-STRIP PCR TUBES



Thermocyclers



Sealing Films



Biologix Group Limited +86-531-67802668

Inquiry@BiologixGroup.com

www.BiologixGroup.com