



**Comments for pBAD/Myc-His A,
4094 nucleotides**

araBAD promoter region: bases 4-276
 Initiation ATG: bases 319-321
 Multiple cloning site: bases 317-370
myc epitope: bases 377-406
 Polyhistidine tag: bases 422-439
rrnB transcription termination region: bases 545-702
 Ampicillin ORF: bases 981-1841
 pBR322 origin: bases 1986-2659
 AraC ORF: bases 4068-3190

* Version A does not contain
SnaB I or *Xba* I.
 Version B contains *Xba* I only.
 Version C contains *SnaB* I only.