

Secure

Advanced

Efficient

# New Members of Cryogenic



### ID-Color™ Cardboard Freezer Boxes

- New hinged lid design eliminates metal rivets and ensures that lids are never misplaced or mixed up
- Alphanumeric labeling featured on the lid, bottom, and four sides enables quick and easy identification of samples
- Sturdy, straight plastic dividers withstand ultra-low temperatures

### 3.75 Inch PC Freezer Boxes

- Constructed of sturdy polycarbonate, for use at temperatures from -196°C to 121°C
- 3.75 inches tall, 81 wells, for 5.0 ml microtubes and cryogenic vials
- 5 assorted colors for color-coding : white, red, yellow, green, and blue

### 2ml Cryogenic Vials

- Suitable for long-term storage in liquid nitrogen, and can be repeatedly frozen and thawed

Ordering Information

#### ID-Color™ Cardboard Freezer Boxes

Cat #	Height	Well	Color	Divider	Unit
90-8281	2 in	81-well	Assorted	Grid Divider	5 Pieces/Pack, 4 Packs/Case
90-8200	2 in	100-well	Assorted	Cylindrical Divider	5 Pieces/Pack, 4 Packs/Case

#### 3.75 Inch PC Freezer Boxes

Cat #	Height	Well	Color	Material	Unit
90-9381	3.75 in	81-well	Assorted	PC	5 Pieces/Pack, 4 Packs/Case

#### 2ml Cryogenic Vials

Cat #	Volume	Bottom	Type	Sterile	Unit
88-321X	2 ml	Self-standing	Internal Thread	Yes	25 Tubes/Bag, 20 Bags/Box, 2 Boxes/Case

\*X - Cap Color: 0 - White, 1 - Red, 2 - Green, 3 - Blue, 4 - Yellow, 5 - Purple