



**Comments for pYES2/CT:  
5963 nucleotides**

*GAL1* promoter: bases 1-451

*GAL1* forward priming site: bases 414-437

T7 promoter/priming site: bases 475-494

Multiple cloning site: bases 501-594

V5 epitope: bases 607-648

Polyhistidine (6xHis) region: bases 658-675

*CYC1* transcription termination signal: bases 708-961

*CYC1* reverse priming site: bases 725-743

pUC origin: bases 1145-1818

Ampicillin resistance gene: bases 1963-2823 (complementary strand)

*URA3* gene: bases 2841-3948 (complementary strand)

2 $\mu$  origin: bases 3952-5423

f1 origin: bases 5491-5946 (complementary strand)