

Eppendorf Certificate

Certificate of Quality

Mastertip®

Mastertip pipette tips are single-use plastic tips which, in combination with the Biomaster® pipette, form a perfectly functioning pipetting system.

Quality Assurance

Mastertip pipette tips are manufactured in a controlled ISO class 8 clean room environment according to ISO 14644-1:2016-06 and are subject to regular quality inspections within a statistical process control.

The Eppendorf Quality Management System is certified according to the international standards ISO 9001:2015-09, ISO 13485:2016-03, and ISO 14001:2015-09.

Materials

Mastertip pipette tips are made of virgin polypropylene and polyethylene of highest purity and quality. Neither the material suppliers of Eppendorf nor Eppendorf use or intentionally incorporate the following agents into the materials Eppendorf is using for the production of Mastertip pipette tips:

- Slip agents (including oleamide, erucamide, stearamide)
- Biocides (including di(2-hydroxyethyl)methyl-dodecylammonium salts (DiHEMDA))
- Plasticizers (softeners/phtalates)

Eppendorf confirms that all plastic materials fulfil the requirements as described in the current version of the FDA, CFR, Title 21 (Food and Drugs) in chapter 178.2010 “Antioxidants and Stabilizers for Polymers” and chapter 177.1520 (a)(1), (b) and (c)1.1 “Olefin polymers” regarding substances used for manufacturing of materials and articles or components of articles intended to come into contact with food.

Mastertip pipette tips and all product components do not contain metallic dyes or latex.

Hamburg, November 2020



Joana Tziolis
Expert Liquid Handling
LH-Consumables



Monika Schneider
Head of Global Quality Management &
Regulatory Affairs

Your local distributor: www.eppendorf.com/contact
Eppendorf AG · Barkhausenweg 1 · 22339 Hamburg · Germany
E-mail: eppendorf@eppendorf.com

ISO 9001
Certified

ISO 13485
Certified

ISO 14001
Certified