




**Product Specification**  
**SARMIX® SM1, Rocking mixer, 230 V, 50 Hz - 60 Hz, for outpatient clinic, doctor's practice and laboratory, with EU plug**



### Product description

Order number	90.180.200
Product description	SARMIX® SM1, Rocking mixer, supply voltage: 230 V, 50 Hz - 60 Hz, for outpatient clinic, doctor's practice and laboratory, with EU plug

### Product characteristics

Design	EU plug
Suitable for	outpatient clinic, doctor's practice and laboratory
Power consumption during operation	4 W
Frequency	50 Hz - 60 Hz

### Size

Width of product	335 mm
Length of product	189 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.





**Product Specification**  
**SARMIX® SM1, Rocking mixer, 230 V, 50 Hz - 60 Hz, for outpatient clinic, doctor's practice and laboratory, with EU plug**

**Purity & certification**

Product in accordance with MDD or IVDD	no medical device   IVD
CE certified	CE
Lot or serial number	no

**Packaging**

Minimum order qty.	1
Piece(s) / inner box	1
Piece(s) / case	1
Piece(s) / pallet	1
Case volume	0,0139 cbm
Weight of product	2,2 kg
Weight of case	2,8 kg
EAN of inner box	4038917101868
EAN case	4038917008549

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

This specification is confidential and the property of Sarstedt. It is neither to be duplicated nor made available to third parties without our prior written consent.

