



Comments for pCR[®]2.1-TOPO[®]
3931 nucleotides

- LacZα* fragment: bases 1-547
- M13 reverse priming site: bases 205-221
- Multiple cloning site: bases 234-357
- T7 promoter/priming site: bases 364-383
- M13 Forward (-20) priming site: bases 391-406
- f1 origin: bases 548-985
- Kanamycin resistance ORF: bases 1319-2113
- Ampicillin resistance ORF: bases 2131-2991
- pUC origin: bases 3136-3809