

Now Available from **BIOLGX**[®]
NEW—Extended 3mL Transfer Pipets!



➤ **More Accurate**

Transparent, thicker pipet walls, evenly distributed thickness, and molded graduation marks ensure consistent pipetting.


➤ **More Reliable**

Constructed of durable, high quality polyethylene material, packaged individually, and gamma radiation sterilized.

➤ **More Convenient**

Integrated pipet and bulb replace awkward rubber bulbs and Pasteur pipets; extended, flexible stems allow easy pipetting from even the deepest solution basins.

What can the extended transfer pipet do?



➤ Replaces rubber bulbs and glass pipets, reducing the likelihood of contamination



➤ Makes pipetting from deeper basins easier



➤ Aspirates and dispenses liquids evenly, without retention

Ordering Information:

Cat #	Package	Length	Graduated To	Sterile	Unit
30-0135	Bulk	160.0mm	1ml	No	500 / Pack, 2000 / Case
30-0135A1	Individually Wrapped	160.0mm	1ml	Yes	500 / Pack, 2000 / Case
30-0138	Bulk	160.0mm	3ml	No	500 / Pack, 2000 / Case
30-0138A1	Individually Wrapped	160.0mm	3ml	Yes	500 / Pack, 2000 / Case
30-0238A1	Individually Wrapped	184.0mm	3ml	Yes	500 / Pack, 2000 / Case



Contact Us: Tel:1-816-6597165 Inquiry@BiologixGroup.com www.BiologixGroup.com