


Retractor frame

Quality - Made in Germany





Picture: Retractor frame with Hot Plate 062

Accessory for:

- Hot Plate 062
- **Hot Plate A3**
- · For variable fixing of objects to magnetic holders
- Stainless steel (magnetic)

Retractor frame for Hot Plate 062

Inner dimensions: 21.0 x 15.0 cm Outer dimensions: 27.7 x 21.7 cm

REF 15943

Retractor frame for Hot Plate A3

Inner dimensions: 34.0 x 22.0 cm Outer dimensions: 42.4 x 30.4 cm

REF 15942

Scope of supply:

- Retractor frame
- Silicone base

(Magnetic holders not in scope of supply)

CE

Rev. 2_02/2020