



**Comments for pTrcHis2/lacZ
7552 nucleotides**

trc promotor region: bases 190-382
 -35 region: bases 193-198
 -10 region: bases 216-221
lac operator (*lacO*): bases 228-248
rrnB antitermination region: bases 264-333
 gene 10 region: bases 346-354
 Ribosome binding site: bases 369-373
 Minicistron ORF: bases 383-409
 Reinitiation RBS: bases 398-403

Expression ATG: bases 413-415
lacZ ORF: bases 467-3523
myc epitope: bases 3617-3649
 Polyhistidine tag: bases 3662-3679
mycHis reverse priming site: bases 3654-3673
 Ampicillin resistance ORF: bases 4222-5082
 pBR322 origin: bases 5227-5900
Lac Repressor (*lacI^q*) ORF: bases 6554-7513