



**Comments for pFLD/CAT:
5014 nucleotides**

FLD1 promoter region: bases 1-597
 CATORF: bases 612-1268
 V5 epitope: bases 1284-1325
 Polyhistidine (6xHis) tag: bases 1335-1352
AOX1 transcription termination region: bases 1363-1691
 3' *AOX1* priming site: bases 1438-1458
bla promoter: bases 2109-2115
 Ampicillin resistance gene: bases 2150-3010
 pUC origin: bases 3155-3827
TEF1 promoter: bases 3828-4230
 T7 promoter: bases 4236-4302
 Zeocin resistance gene: bases 4303-4677
CYC1 transcription termination region: bases 4678-4995