



**Comments for pEXP2-GW/lacZ
5865 nucleotides**

T7 promoter: bases 20-36
 T7 promoter/priming site: bases 20-39
attB1 site: bases 98-122
lacZ gene: bases 164-3220
attB2 site: bases 3240-3264
 V5 epitope: bases 3284-3325
 V5 (C-term) reverse priming site: bases 3293-3313
 Polyhistidine (6xHis) region: bases 3335-3352
 T7 transcription termination region: bases 3390-3474
 Zeocin™ resistance gene: bases 3618-3992
 Ampicillin resistance gene: bases 4013-4873
 pUC origin: bases 5018-5691