

BIOL@GLX[®]

Chambered Cell Culture Slides

Cat#	Well	Material Chamber/ Slide/ Holder	Color	Growth Area /Well(cm ²)	Working Vol. /Well(ml)	Sterile	Pakaging
07-2101	1	PS /Glass /PP	Clear	9.40	2.50-5.50	Yes	6 /Pack , 12 /Case
07-2102	2	PS /Glass /PP	Clear	4.55	1.20-2.50	Yes	6 /Pack , 12 /Case
07-2104	4	PS /Glass /PP	Clear	2.13	0.50-1.30	Yes	6 /Pack , 12 /Case
07-2108	8	PS /Glass /PP	Clear	0.98	0.20-0.60	Yes	6 /Pack , 12 /Case

Don't forget Biologix's other Cell Culture Products



Contact Us: Tel:1-816-6597165

Inquiry@BiologixGroup.com

www.BiologixGroup.com