



Comments for:

CMV promoter:
 T7 promoter/priming site:
attB1 site:
 CAT gene:
attB2 site:
 EmGFP/YFP:
 TK polyA reverse priming site:
 TK polyadenylation signal:
f1 origin:
 SV40 early promoter and origin:
 EM7 promoter:
 Blasticidin resistance gene:
 SV40 early polyadenylation signal:
 pUC origin (c):
 Ampicillin (*bla*) resistance gene (c):
bla promoter (c):

(c) = complementary strand

**pcDNA™ 6.2/
 C-EmGFP/GW/CAT
 6448 nucleotides**

**cDNA™ 6.2/
 C-YFP/GW/CAT
 6448 nucleotides**

102-689
 733-752
 785-809
 813-1472
 1471-1495
 1509-2228
 2250-2268
 2243-2514
 2550-3009
 3005-3313
 3368-3434
 3435-3833
 3991-4121
 4504-5177
 5322-6182
 6177-6281

102-689
 733-752
 785-809
 813-1472
 1471-1495
 1509-2228
 2250-2268
 2243-2514
 2550-3009
 3005-3313
 3368-3434
 3435-3833
 3991-4121
 4504-5177
 5322-6182
 6177-6281