


**Comments for pTracer™-CMV
5717 nucleotides**



CMV promoter: bases 209-863
T7 promoter priming site: bases 863-882
Multiple cloning site: bases 901-1010
pcDNA3.1 reverse priming site: bases 1022-1039
BGH polyadenylation signal: bases 1021-1235
f1 origin: bases 1298-1711
SV40 promoter and origin: bases 1776-2101
EM-7 promoter: bases 2117-2183
GFP-Zeocin™ resistance gene: bases 2184-3260
 GFP ORF: bases 2184-2888
 Zeocin™ ORF (no Met): bases 2889-3260
SV40 polyadenylation: bases 3269-3471
ColE1 origin: bases 3903-4576
Ampicillin resistance gene: bases 4721-5580



The sequence of pTracer™-CMV has been compiled from information in sequence databases, published sequences, and other sources. This vector has not been completely sequenced. If you suspect an error in the sequence, please contact Invitrogen's Technical Service Department.