

BRC-5180UT High Speed Bench-Top Touch Panel Centrifuge Rotor Types



Order No.	Rotor Type	Max. Speed (rpm)	Max. RCF (xg)	Max. Capacity
INO518UT-1	Swing Rotor (Round Bucket)	4000	2580	4x180ml
INO518UT-2	Microplate Rotor (Fame Type)	4000	2020	2x3x96 well
INO518UT-3	Microplate Rotor (Box Type)	4000	1980	2x3x96 well
INO518UT-4	Swing Rotor	5000	4420	4x1x50ml
INO518UT-5	Swing Rotor	5000	4420	4x1x100ml
INO518UT-6	Swing Rotor	4000	2820	4x2x50ml
INO518UT-7	Swing Rotor	4000	2820	4x4x15ml
INO518UT-8	Swing Rotor	4000	2580	4x6x10/7ml
INO518UT-9	Swing Rotor	4000	2170	4x6x5ml
INO518UT-10	Swing Rotor	4000	2580	4x4x10/7ml
INO518UT-11	Swing Rotor	4000	2170	4x4x5ml
INO518UT-12	Fixed Rotor	16000	15760	4x8 PCR
INO518UT-13	Fixed Rotor	15000	21420	6x8 PCR
INO518UT-14	Fixed Rotor	16000	17480	8x8 PCR
INO518UT-15	Fixed Rotor	15000	22930	12x8 PCR
INO518UT-16	Fixed Rotor	15000	22920	40x0.5ml
INO518UT-17	Fixed Rotor	21000	30910	12x1.5ml
INO518UT-18	Fixed Rotor	16000	23440	24x1.5ml
INO518UT-19	Fixed Rotor	14000	20800	30x1.5ml
INO518UT-20	Fixed Rotor	13000	17930	48x1.5ml
INO518UT-21	Fixed Rotor	16000	22020	16x5ml
INO518UT-22	Fixed Rotor	16000	21500	6x10ml
INO518UT-23	Fixed Rotor	15000	22680	12x10ml
INO518UT-24	Fixed Rotor	13000	19490	16x10ml
INO518UT-25	Fixed Rotor	13000	17790	8x15ml (Conical)
INO518UT-26	Fixed Rotor	11000	14330	12x15ml
INO518UT-27	Fixed Rotor	5000	3500	24x15ml
INO518UT-28	Fixed Rotor	5000	3830	30x15ml
INO518UT-29	Fixed Rotor	14000	19060	6x30ml
INO518UT-30	Fixed Rotor	13000	18840	6x50ml (Conical)
INO518UT-31	Fixed Rotor	13000	18730	6x50ml (Round)
INO518UT-32	Fixed Rotor	5000	3860	12x50ml
INO518UT-33	Fixed Rotor	13000	18940	4x85ml
INO518UT-34	Fixed Rotor	12000	14850	4x100ml
INO518UT-35	Fixed Rotor	10000	11380	6x100ml
INO518UT-36	Fixed Rotor	12000	15800	24 Capillaries