

BRC-5180T High Speed Bench-Top Touch Panel Centrifuge Rotor Types



Order No.	Rotor Type	Max. Speed (rpm)	Max. RCF (xg)	Max. Capacity
INO518T-1	Swing Rotor (Round Bucket)	4000	2580	4x180ml
INO518T-2	Microplate Rotor (Fame Type)	4000	2020	2x3x96 well
INO518T-3	Microplate Rotor (Box Type)	4000	1980	2x3x96 well
INO518T-4	Swing Rotor	5000	4420	4x1x50ml
INO518T-5	Swing Rotor	5000	4420	4x1x100ml
INO518T-6	Swing Rotor	4000	2820	4x2x50ml
INO518T-7	Swing Rotor	4000	2820	4x4x15ml
INO518T-8	Swing Rotor	4000	2580	4x6x10/7ml
INO518T-9	Fixed Rotor	14000	12070	4x8 PCR
INO518T-10	Fixed Rotor	13000	16080	6x8 PCR
INO518T-11	Fixed Rotor	14000	13390	8x8 PCR
INO518T-12	Fixed Rotor	13000	17220	12x8 PCR
INO518T-13	Fixed Rotor	14000	19960	40x0.5ml
INO518T-14	Fixed Rotor	18500	23990	12x1.5ml
INO518T-15	Fixed Rotor	15000	20600	24x1.5ml
INO518T-16	Fixed Rotor	13500	19340	30x1.5ml
INO518T-17	Fixed Rotor	11000	12840	48x1.5ml
INO518T-18	Fixed Rotor	15000	19350	16x5ml
INO518T-19	Fixed Rotor	14000	16460	6x10ml
INO518T-20	Fixed Rotor	13000	17030	12x10ml
INO518T-21	Fixed Rotor	11000	13950	16x10ml
INO518T-22	Fixed Rotor	11000	12730	8x15ml
INO518T-23	Fixed Rotor	9000	9590	12x15ml
INO518T-24	Fixed Rotor	5000	3500	24x15ml
INO518T-25	Fixed Rotor	5000	3830	30x15ml
INO518T-26	Fixed Rotor	12000	14000	6x30ml
INO518T-27	Fixed Rotor	12000	16050	6x50ml (Conical)
INO518T-28	Fixed Rotor	12000	15960	6x50ml (Round)
INO518T-29	Fixed Rotor	11000	13560	4x85ml
INO518T-30	Fixed Rotor	10000	10310	4x100ml
INO518T-31	Fixed Rotor	8000	7280	6x100ml
INO518T-32	Fixed Rotor	12000	15800	24 Capillaries
INO518T-33	Fixed Rotor	15000	19920	4x5ml
INO518T-34	Fixed Rotor	16000	16450	16x5ml
INO518T-35	Fixed Rotor	14000	19750	8x30ml