

INOIA-50 Quantitative Immunoassay Analyzer



Easy to use

Display : 5.0" touch screen
Calibration Mode : Magnetic card
Maintenance-free

Rapid Detection

Fast mode and standard mode
Rapid Detection : Results available within 3-15min
Measurement Time : 10sec

Technical Specifications

Principle : Immunofluorescence
Sample Type : Whole blood, serum, plasma
Sample Volume : 100ul
Storage : 50,000 results storage
Output : Internal thermal printer
Support LIS and Barcode scanner
Update : 2 USB ports
Work Temperature : 10 - 30°C
Power Supply : AC 100/240V, 50/60Hz
Net Weight : 4 Kg



Magnetic Card Calibration



Pipette Sample



Insert Cassette



Result

Test Menu

Item	Applications	Item	Applications
PCT	Inflammation	T3	Thyroid
hs-CRP	Inflammation	T4	Thyroid
IL-6	Inflammation	TSH	Thyroid
cTnI	Myocardial Infarction	LH	Fertility
Myo	Myocardial Infarction	FSH	Fertility
CK-MB	Myocardial Infarction	PRL	Fertility
cTnI/CK-MB/Myo	Myocardial Infarction	PROG	Fertility
D-Dimer	Coagulation	β -HCG	Fertility
HbA1c	Diabetes	TESTO	Fertility
NT-proBNP	Heart Failure		